

Product

PrepMan[®] Ultra (50mL)
(Sample Preparation Reagent)

Applied Biosystems Product Number

4318930

0.5 mL Microtubes with Caps
(500/bag)

VWR Part Number

16466-050

Preparation

1. Aliquot 100 µL PrepMan[®] Ultra in 1.5mL or 0.5 mL microtubes with caps
2. Suspend a small loop of cell mass in microtube and vortex for 15 seconds
3. Heat sample at 95°C for a minimum of 10 minutes in a water bath or heat block
4. Cool, enclose in secondary container (zip closure plastic bag) and ship at ambient temperature to Accugenix, Inc.:

223 Lake Drive
Newark, DE 19702-9912
1.800.886.9654
www.accugenix.com

5. Complete the "Description of Goods" section on the shipper's paperwork (Airway Bill and Commercial Invoice) with the description:

"non-hazardous, sterile DNA suspensions not restricted by IATA/ IMDG/ DOT"