

Customer:	Accugenix, Inc.	Date of Report:	2010-03-11
Cust Address:	223 Lake Drive, Newark, DE, 19702, USA	Cust Account#:	111555 (ACC2)
Cust Sample ID:	S.epidermidis ATCC 14990	Accugenix C#:	C1000116
Orig. Due Date:	2009-11-09	IRF#	100012

Identification: **Staphylococcus epidermidis**
Confidence Level: **Species**

Quality	Top Species Match(es) for Organism	Score	Run Date/Time
+++	Staphylococcus epidermidis	2.336	2009-11-03 11:47:57

Rank (Quality)	Top Ten Matched Patterns	Score	NCBI ID
1 (+++)	Staphylococcus epidermidis ATCC 14990 THL	2.336	1282
2 (+++)	Staphylococcus epidermidis ATCC 14990 THL	2.336	1282
3 (+++)	Staphylococcus epidermidis ATCC 12228 THL	2.077	1282
4 (+++)	Staphylococcus epidermidis ATCC 12228 THL	2.077	1282
5 (+++)	Staphylococcus epidermidis 10547 CHB	2.074	1282
6 (+++)	Staphylococcus epidermidis 10547 CHB	2.074	1282
7 (+++)	Staphylococcus epidermidis DSM 3269	1.967	1282
8 (+++)	Staphylococcus epidermidis DSM 3269	1.967	1282
9 (+++)	Staphylococcus epidermidis ATCC 12228 CHB	1.860	1282
10 (+++)	Staphylococcus epidermidis ATCC 12228 CHB	1.860	1282

Meaning of Score Values

Range	Description	Symbols	Color
1.750 ... 3.000	probable species identification	(+++)	green
0.000 ... 1.749	not a reliable species identification	(-)	red

Prepared by: E 11 Mar 10
 Reviewed by: SS 11 Mar 10
 QA Approved by: M 11 MAR 10